

Laundry Instructions

Vinyl Holsters & Mattress Blocks

Your vinyl holsters (including Bedpan Holster, Fracture Pan Holster, Urinal Holster, Wheelchair Urinary Drainage Bag Holster and the Tube Hider) can be machine washed in warm water. Tumble dry on low and remove promptly. **DO NOT USE BLEACH.**

The Mattress Block can be sanitized by wiping, just like a mattress. **DO NOT USE BLEACH.**

Questions?

Problems?

Suggestions?

Call 1-800-336-1414

Powderhorn Industries, 931 North Park Ave., P.O. Box 1443, Montrose, CO 81402